# MASCOT Distiller : 機能

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Mascot Distiller - QEP1\_SpikeIn\_230914\_1\_3ng\_1data\*

File Edit View Processing **Analysis** Tools Windows Help

Mascot Search All Peaklists...

Denovo Search Current peaklist...  
Digest Protein... Mascot Search (All Tags)...  
Fragment Peptide... Current Tag...  
Analysis Info All Tags in Current Peaklist...  
Calculate XIC...

Peak Lists

Grouped Scans  
MS/MS 12783, m/z 483.452, 5+ 5920:  
Grouped Scans  
MS/MS 12787, m/z 432.732, 2+ 5921:  
Grouped Scans  
MS/MS 12788, m/z 864.389, 2+ 5922:  
Grouped Scans  
MS/MS 12789, m/z 366.725, 2+ 5923:  
Grouped Scans  
MS/MS 12790, m/z 366.725, 2+ 5923:

[6] Analysis  
-> Mascot Search  
-> All Peaklists

Mascot Search (localhost)

### MASCOT MS/MS Ions Search

Your name: takaesu Email: [ ]

Search title: C:\temp\replicatesDataset\QEP1\_SpikeIn\_230914\_1\_3ng\_270914.raw

Database(s): SwissProt (AA)

Amino acid (AA): AmbystomaMexicanum, Caudata, contaminants, cRAP, IPI\_human, M\_smeigmatis\_2, M\_smeigmatis\_ATCC700084, Mouse, NCBIprot

Taxonomy: All entries

Enzyme: Trypsin Allow up to: 1 missed cleavages

Quantitation: None

Crosslinking: None

Fixed modifications: Carbamidomethyl (C)

Variable modifications: Oxidation (M)

Display all modifications:

Peptide tol. ±: 0.1 Da # <sup>13</sup>C: 0 MS/MS tol. ±: 0.1 Da

Peptide charge: 2+ Monoisotopic:  Average

Data file: [ ]

Data format: Mascot generic Precursor: [ ] m/z

Instrument: ESI-TRAP Error tolerant:

Decoy:



Start Search ... Reset Form

Copyright © 2018 Matrix Science Ltd. All Rights Reserved.

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Status | Event Log | Task Editor | **Parameter Editor**

Parameter set

Filename: C:\ProgramData\Matrix Science\Mas...default.par [New] [Open ...] [Save] [Save As ...]

All Searches

User name <mascot\_user\_full\_name> User email <mascot\_user\_email>

Search title <taskname> (<parameters>), submitted from Daemon on <localhost>

Taxonomy ..... Homo sapiens (human) Report top AUTO hits

Databases SwissProt [Select Databases] Protein mass kDa

Decoy  Enzyme Trypsin

Fixed modifications Carbamidomethyl (C) [Select Modifications] Max. missed cleavages 1

Variable modifications Oxidation (M) Monoisotopic  Average  Peptide charge 2+ and 3+

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# Daemon + Distiller の操作手順 : Task Editor

The screenshot displays the Mascot Daemon Task Editor window. The interface includes a menu bar (File, Edit, Help), a status bar (Status, Event Log, Task Editor, Parameter Editor), and several main sections: Owner, Task (with a text input field containing 'search test 01' and 'Run' and 'New' buttons), Parameter set (with a file path and a selection button), Data file list (with a list of files and an 'Add Files ...' button), Data import filter (with a dropdown menu set to 'Mascot Distiller' and an 'Options ...' button), Schedule (with radio buttons for 'Start now', 'Start at', 'Start on completion of', 'Real-time monitor', and 'Follow-up'), Actions (with 'Auto-export ...' and 'External processes ...' buttons), and Follow-up (with a dropdown menu set to 'No follow-up required' and a '0' input field).

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The screenshot shows the 'Mascot Distiller data import options' dialog box. The 'Data File Format' dropdown is set to 'Thermo Xcalibur'. The 'Mascot Distiller Processing Options' field contains the path 'C:\ProgramData\Matrix S...%prof\_prof.ThermoXcalibur.opt'. The 'Output MS/MS Fragments as' section has 'MH+' selected. The 'Distiller Project' section has 'Save' checked. The 'Scan level parameters' section has 'Output INSTRUMENT=' checked and set to 'Default' for 'Unknown' scans.

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Re-search  All  Non-significant  Unassigned [\[help\]](#) Export As

Not what you expected? Try [the select summary](#).

- Search parameters
- Score distribution
- Modification statistics for all protein families
- Legend

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