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The screenshot shows a web browser window with the URL <https://www.matrixscience.com/licensing/register?sid=DLLIC>. The page title is "Mascot Licence Registration".

Instructions at the top of the form state: "Please type in your contact details to finish creating the licence." and provide a list of notes: "If any of the product and licence details shown below are incorrect, please contact Matrix Science support." and "Required fields are marked with an asterisk\*." It also notes that the licence file will be sent to the email address provided and advises adding support@matrixscience.com to a spam filter's whitelist.

The form fields are as follows:

- Product key:** [Redacted]
- Product family:** Mascot Server
- Licence type:** Permanent
- Licensed features:** Version 3.0, 1 CPUs, 6
- Email address\*:** takaesu@matrixscience.com
- Verify email address\*:** takaesu@matrixscience.com
- Full name\*:** Hirotomo Takaesu
- Organisation\*:** Matrix Science Japan
- Address line 1\*:** 1-6-10,Higasi-Ueno
- Address line 2:** [Empty]
- City\*:** Tokyo
- Zip/postcode:** 1100015
- State/region:** [Empty]
- Country\*:** Japan
- Phone number:** [Empty]
- Fax number:** [Empty]

Below the form, there is a section for "Email address (CC)" with two radio button options:

- Matrix Science may send me technical support information by email, such as the availability of new releases, upgrades, and support contracts.
- I do not wish to receive any communications from Matrix Science.

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**MASCOT PROTEIN IDENTIFICATION SYSTEM**

**End-user Licence Agreement**

**IMPORTANT – PLEASE READ CAREFULLY:** This End User Licence Agreement is a legally binding contract between you (either an individual or a single corporate entity) and Matrix Science Limited for the product identified above, which includes computer software, electronic documentation, any printed documentation, and any subsequent updates and supplements (the “Software”).

By installing or using the Software, you agree to be bound by the terms of this agreement. If you do not agree to the terms of this agreement, we are unwilling to license the Software to you. In this case, do not install or use the Software. Return the package that included the Software to Matrix Science Limited or their authorised distributor within 30 days of receipt for a full refund.

**1 Licence**

Matrix Science Limited owns the copyright in the Software and all other copies which you receive with the package and all other copies which you receive with the package.

This licence is personal to you (either an individual or a single corporate entity) as the purchaser of a licence to use the Software and the licence granted herein is for your benefit only.

You may not use the Software in any way that permits unlicensed access to the Software. In particular, individuals who are not party to this licence or the general

I do not accept the terms of the licence agreement.



I accept the terms of the licence agreement.\*

Create Licence

Copyright © 1998-2024 Matrix Science Ltd. All Rights Reserv

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Name = PRIDE\_Contaminants\_ Family = /usr/local/mas  
Filename = PRIDE\_Contaminants\_20160908.msp Pathname = /usr  
Status = In use [Statistics](#) [Compi](#)  
State Time = Tue Oct 22 16:02:43 # searches = 0  
Mem mapped = NO Request to mem map = NO Request unmap = NO  
Number of threads = 1 Current = YES Type = Spectral Libra

Name = SwissProt\_ Family = /usr/local/mas  
Filename = SwissProt\_2024\_05.fasta Pathname = /usr/local  
Status = In use [Statistics](#) [Unid](#)  
State Time = Tue Oct 22 16:01:32 # searches = 0  
Mem mapped = YES Request to mem map = YES Request unmap =  
Number of threads = -1 Current = YES Type = Amino acid

Licence information

localhost/mascot/x-cgi/ms-status.exe?Show=VIEWLIC

### Mascot Server Licence Information

[Register a new product key](#) [View database status](#) [Reload this page](#)

Please include *all* the contents of this page when requested to provide this information to technical support.

Mascot Server version: 3.0.0  
Licence path: /usr/local/mascot/config/Licdb


Licence(s) found:

Product Key	Start	End	Status
Active:			
	-10-27 *		OK
Feature: Mascot Server - Core functionality (v3.0)			
Feature: Mascot Server - CPU units (1)			
Feature: Mascot Server - Cores per CPU units (6)			
Company: Matrix Science Japan			
User: [REDACTED]			
Distributor: Matrix Science KK			
Inactive:			

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
Latest predefined definitions files are from Thu Mar 8 20:27:25 2018 (FASTA databases: databases\_1.xml) and Tue Nov 8 00:24:22 2016 (spectral libraries: libraries\_1.xml).  
Full database status is available on [the database status page](#).

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The screenshot shows the 'Databases and spectral libraries' page in the Database Manager. The 'Settings' link in the left sidebar is highlighted with a red box. A yellow arrow points from this link to the 'Proxy settings' section below. In the 'Proxy settings' section, the test results for HTTP and HTTPS are both 'succeeds', which are also highlighted with red boxes. A yellow box contains the text: '现阶段でインターネット接続できるかのテスト結果。接続できていれば「Succeeds」'. The 'Edit proxy settings' button is also highlighted with a red box, and a yellow box below it says: '設定画面は「Edit proxy settings」クリック'. A large yellow arrow at the bottom of the page points downwards.

**Database Manager**

Databases (2)

Parse rules (7)

Scheduled updates (0)

Running tasks (0)

**Settings**

fasta

able predefined  
nition

chronise custom  
nitions

ate new

rary

### Databases and spectral libraries

Name	Mode ?	Type ?	Status	Latest task
PRIDE_Contaminants	predefined	SL	In use	Update succeeded (view log)
SwissProt	predefined	AA	In use	Update succeeded (view log)

Get new files Deactivate

Latest predefined definitions files are from Mon Jan 29 22:16:40 2024 (FASTA databases: databases\_1.xml) and Thu Oct 13 21:39:16 2022 (spectral libraries: libraries\_1.xml).

Full database status is available on the [database status page](#).

Refresh

Disable

### Proxy settings

You may need to define proxy settings to enable Database Manager to access the Internet, depending on your local network and Internet connection.

Proxy settings defined in:

- HTTP proxy strate
- HTTPS proxy stra
- FTP proxy strateg

HTTP test **succeeds**; can access [http://http-connectivity-test.matrixscience.com/db\\_manager.txt](http://http-connectivity-test.matrixscience.com/db_manager.txt) from the Mascot server machine.

HTTPS test **succeeds**; can access [https://https-connectivity-test.matrixscience.com/db\\_manager.txt](https://https-connectivity-test.matrixscience.com/db_manager.txt) from the Mascot server machine.

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
The screenshot shows the 'Proxy settings' section of the Mascot Database Manager. The page title is 'Mascot Database Manager Proxy settings'. Under 'mascot.dat settings', it states: 'The recommended proxy type is Auto. Try in order, until one succeeds.' The 'Proxy type' dropdown menu is open, showing options: Auto, None, Environment, Registry, Specify, and WPAD. The 'Auto' option is highlighted in blue. A yellow callout box points to this menu with the text: 'この画面で proxy 設定する場合、“Specify”を選択'. Below the dropdown, the 'Proxy type' is set to 'Auto'. The 'Host (HTTP proxy) (?)' field is empty, with a yellow callout box pointing to it: '← Proxy サーバアドレス 例: http://proxy.mascot.jp/proxy'. The 'Port (HTTP proxy) (?)' field is empty, with a yellow callout box: '← Proxy ポート。例: 8080'. The 'User (HTTP proxy) (?)' and 'Password (HTTP proxy) (?)' fields are empty, with a yellow callout box: '← 認証ユーザー名とパスワード。'. Below these fields, there are two checkboxes: 'Use HTTP proxy for FTP (?)' and 'Use HTTP proxy for HTTPS (?)', both of which are unchecked. A red box highlights these checkboxes. Below the checkboxes is a 'Save' button, which is also highlighted with a red box. A yellow callout box points to the 'Save' button: '← 設定完了後、Save'. At the bottom of the page, there is a 'Back' button. The page also contains test results: 'HTTP test succeeds; can access http://http-connectivity-test.matrixscience.com/db\_manager.txt from the Mascot server.' and 'HTTPS test succeeds; can access https://https-connectivity-test.matrixscience.com/db\_manager.txt from the Mascot server.' (FTP proxy settings not tested.).

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[https://www.matrixscience.com/help/seq\\_db\\_setup.html](https://www.matrixscience.com/help/seq_db_setup.html)

[https://www.matrixscience.com/help/seq\\_db\\_setup\\_library.html](https://www.matrixscience.com/help/seq_db_setup_library.html)

Name	
contaminants	<input type="button" value="Enable"/>
cRAP	<input type="button" value="Enable"/>
Environmental_EST	<input type="button" value="Enable"/>
Fungi_EST	<input type="button" value="Enable"/>
Human_EST	<input type="button" value="Enable"/>
Invertebrates_EST	<input type="button" value="Enable"/>
Mammals_EST	<input type="button" value="Enable"/>
Mus_EST	<input type="button" value="Enable"/>
NCBIprot	<input type="button" value="Enable"/>
neXtProt	<input type="button" value="Enable"/>
Plants_EST	<input type="button" value="Enable"/>
Prokaryotes_EST	<input type="button" value="Enable"/>
Rodents_EST	<input type="button" value="Enable"/>
SwissProt_AC	<input type="button" value="Enable"/>
SwissProt_ID	<input type="button" value="Enable"/> Already
Trembl_AC	<input type="button" value="Enable"/>
Trembl_ID	<input type="button" value="Enable"/>
UniRef100	<input type="button" value="Enable"/>
UP186698_X_jaevis	<input type="button" value="Enable"/>
UP1940_C_elegans	<input type="button" value="Enable"/>
UP2195_D_discoideum	<input type="button" value="Enable"/>
UP2311_S_cerevisiae	<input type="button" value="Enable"/>
UP2485_S_pombe	<input type="button" value="Enable"/>
UP2494_R_norvegicus	<input type="button" value="Enable"/>
UP437_D_rerio	<input type="button" value="Enable"/>
UP5226_T_rubripes	<input type="button" value="Enable"/>
UP5640_H_sapiens	<input type="button" value="Enable"/>
UP589_M_musculus	<input type="button" value="Enable"/>
UP59680_O_sativa	<input type="button" value="Enable"/>
UP625_E_coli_K12	<input type="button" value="Enable"/>
UP6548_A_thaliana	<input type="button" value="Enable"/>
UP6906_C_reinhardtii	<input type="button" value="Enable"/>
UP7305_Z_mays	<input type="button" value="Enable"/>
UP803_D_melanogaster	<input type="button" value="Enable"/>
UP808_M_pneumoniae	<input type="button" value="Enable"/>
UP9136_B_taurus	<input type="button" value="Enable"/>
Vertebrates_EST	<input type="button" value="Enable"/>

Name	
NIST_BSA_IonTrap	<input type="button" value="Enable"/>
NIST_C.elegans_IonTrap	<input type="button" value="Enable"/>
NIST_Chicken_IonTrap	<input type="button" value="Enable"/>
NIST_D.rerio_IonTrap	<input type="button" value="Enable"/>
NIST_Drosophila_IonTrap	<input type="button" value="Enable"/>
NIST_E.coli_IonTrap	<input type="button" value="Enable"/>
NIST_HSA_IonTrap	<input type="button" value="Enable"/>
NIST_Human_HCD	<input type="button" value="Enable"/>
NIST_Human_HCD_2_good	<input type="button" value="Enable"/>
NIST_Human_HCD_3_semitryp	<input type="button" value="Enable"/>
NIST_Human_HCD_ITRAQ_1	<input type="button" value="Enable"/>
NIST_Human_HCD_ITRAQ_2	<input type="button" value="Enable"/>
NIST_Human_HCD_ITRAQ_Phospho	<input type="button" value="Enable"/>
NIST_Human_IonTrap	<input type="button" value="Enable"/>
NIST_Mouse_HCD	<input type="button" value="Enable"/>
NIST_Mouse_HCD_ITRAQ	<input type="button" value="Enable"/>
NIST_Mouse_HCD_ITRAQ_Phospho	<input type="button" value="Enable"/>
NIST_Mouse_IonTrap	<input type="button" value="Enable"/>
NIST_Rat_IonTrap	<input type="button" value="Enable"/>
NIST_Rat_QToF	<input type="button" value="Enable"/>
NIST_S.cerevisiae_IonTrap	<input type="button" value="Enable"/>
NIST_S.cerevisiae_QToF	<input type="button" value="Enable"/>
NIST_S.pombe_IonTrap	<input type="button" value="Enable"/>
NIST_UPS1_IonTrap	<input type="button" value="Enable"/>
PRIDE_Arabidopsis	<input type="button" value="Enable"/>
PRIDE_B.burgdorferi	<input type="button" value="Enable"/>
PRIDE_B.cereus	<input type="button" value="Enable"/>
PRIDE_B.subtilis	<input type="button" value="Enable"/>
PRIDE_C.elegans	<input type="button" value="Enable"/>
PRIDE_Contaminants	<input type="button" value="Enable"/> Already set t
PRIDE_Drosophila	<input type="button" value="Enable"/>
PRIDE_E.coli	<input type="button" value="Enable"/>
PRIDE_G.sulfurreducens	<input type="button" value="Enable"/>
PRIDE_Human	<input type="button" value="Enable"/>
PRIDE_Mouse	<input type="button" value="Enable"/>
PRIDE_Rat	<input type="button" value="Enable"/>
PRIDE_S.cerevisiae	<input type="button" value="Enable"/>
PRIDE_S.typhimurium	<input type="button" value="Enable"/>
PRIDE_Soybean	<input type="button" value="Enable"/>

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**Enable predefined definition**

Library name: NIST\_Human\_IonTrap

Base directory (?): /usr/local/mascot/sequence

Library files will be located in the subdirectory NIST\_Human\_IonTrap of the base directory. If the directory does not already exist, it will be created.

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**Enable predefined definition**

Library name: NIST\_Human\_IonTrap

Sequence directory: /usr/local/mascot/sequence

Reference database: SwissProt

Please choose a reference data source (FASTA file) from the list below. The reference data source is the one displayed for a library hit.

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Database type: Spectral library

Database directory: /usr/local/mascot/sequence/NIST\_Human\_IonTrap/current

Filename pattern: NIST\_Human\_IonTrap\_\*.msp

About to download or copy files

Database files need to be downloaded or copied before database configuration can be completed.

Download task is in queue position 1; waiting for downloading to start. See the task queue for more details. (Download task added Mon Oct 28 16:48:47 2024.)

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The screenshot shows the Mascot Database Manager interface. The left sidebar has a menu with 'Running tasks (1)' highlighted in a red box. An orange arrow points from this menu item to the 'Database Manager tasks' section. In this section, a task is shown with a progress bar at 6.0%. A yellow box with the text '実行待ち→' points to the 'Task process control' section. Another yellow box with '実行中 →' points to the task row. A third yellow box with '実行終了→' points to the 'Task log' section. In the 'Task log', a task is shown with a 'Last message' column containing the text '(success) 'NIST\_Human\_IonTrap' successfully updated.' This text is highlighted with a red box. A large yellow box with the text '(Success)' DB名' successfully updated と表示されたら成功' is overlaid on the 'Last message' column.

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The screenshot shows the Mascot Database Status page with the following entries:

Name	Family	Filename	Status
PRIDE Contaminants	C:/inetpub/mascot/sequence/PRIDE_Contaminants/current/PRIDE_Contaminants/current/PRIDE_Contaminants_20160906.msp	PRIDE_Contaminants_20160906.msp	Not in use
PRIDE Contaminants	C:/inetpub/mascot/sequence/PRIDE_Contaminants/current/PRIDE_Contaminants/current/PRIDE_Contaminants_20160908.msp	PRIDE_Contaminants_20160908.msp	In use
SwissProt	C:/inetpub/mascot/sequence/SwissProt/current/SwissProt_*.fasta	SwissProt_2016_10.fasta	Not in use
SwissProt	C:/inetpub/mascot/sequence/SwissProt/current/SwissProt_*.fasta	SwissProt_2017_01.fasta	In use

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