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- ☒ Web Server
 - ☒ Common HTTP Features
 - ☒ Default Document
 - ☒ Directory Browsing
 - ☒ HTTP Errors
 - ☒ Static Content
 - ☐ HTTP Redirection
 - ☐ WebDAV Publishing
 - ☒ Health and Diagnostics
 - ☒ HTTP Logging
 - ☐ Custom Logging
 - ☐ Logging Tools
 - ☐ ODBC Logging
 - ☐ Request Monitor
 - ☐ Tracing
 - ☒ Performance
 - ☒ Static Content Compression
 - ☐ Dynamic Content Compression
 - ☒ Security
 - ☒ Request Filtering
 - ☐ Basic Authentication
 - ☐ Centralized SSL Certificate Support
 - ☐ Client Certificate Mapping Authentication
 - ☐ Digest Authentication
 - ☐ IIS Client Certificate Mapping Authentication
 - ☐ IP and Domain Restrictions
 - ☐ URL Authorization
 - ☐ Windows Authentication
 - ☒ Application Development
 - ☐ .NET Extensibility 3.5
 - ☐ .NET Extensibility 4.5
 - ☐ Application Initialization
 - ☐ ASP
 - ☐ ASP.NET 3.5
 - ☐ ASP.NET 4.5
 - ☒ CGI
 - ☐ ISAPI Extensions
 - ☐ ISAPI Filters
 - ☐ Server Side Includes
 - ☐ WebSocket Protocol
 - ☐ FTP Server
 - ☐ IIS Hostable Web Core
 - ☒ Management Tools
 - ☒ IIS Management Console
 - ☒ IIS 6 Management Compatibility
 - ☒ IIS 6 Metabase Compatibility
 - ☒ IIS 6 Management Console
 - ☒ IIS 6 Scripting Tools
 - ☒ IIS 6 WMI Compatibility
 - ☒ IIS Management Scripts and Tools
 - ☐ Management Service

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The screenshot shows the 'Mascot Licence Registration' web form. The browser address bar shows 'www.matrixscience.com'. The form includes a search bar and navigation links. A message box states: 'Please type in your contact details to finish creating the licence. If any of the product and licence details shown below are incorrect, please contact Matrix Science support. Required fields are marked with an asterisk*. The licence file will be sent to the email address you supply below. Make sure your spam filter does not reject this message by adding support@matrixscience.com to your [whitelist] or address book.'

The form fields are as follows:

- Product key**: [Text box]
- Product family**: Mascot Server
- Licence type**: Permanent
- Licensed features**: Version [Text box], Version [Text box], 2 CPUs
- Old product key**: [Text box]
- Email address***: takaesu@matrixscience.com
- Verify email address***: takaesu@matrix
- Full name***: Hiroto Tom Takae
- Organisation***: Matrix Science Japan
- Address line 1***: 1-6-10, Higashi-Ueno
- Address line 2**: [Text box]
- City***: Tokyo
- Zip/postcode**: 1100015
- State/region**: Tokyo
- Country***: Japan
- Phone number**: [Text box]
- Fax number**: [Text box]

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Note: In either case, Matrix Science will not use your contact information to send unrelated advertising, and will not disclose your details to any third party.

MASCOT PROTEIN IDENTIFICATION SYSTEM

End-user Licence Agreement

IMPORTANT – PLEASE READ CAREFULLY: This End User Licence Agreement is a legally binding contract between you (either an individual or a single corporate entity) and Matrix Science Limited for the product identified above, which includes computer software, electronic documentation, any printed documentation, and any subsequent updates and supplements (the "Software").

By installing or using the Software, you agree to be bound by the terms of this agreement. If you do not agree to the terms of this agreement, we are unwilling to license the Software to you. In this case, do not install or use the Software. Return the package that included the Software to Matrix Science Limited or their authorised distributor within 30 days of receipt for a full refund.

1 Licence

Matrix Science Limited owns the copyright in the Software contained within this package and all other copies which you are authorised by this agreement to make.

This licence is personal to you (either an individual or a single corporate entity) as the purchaser of a

☐ I do not accept the terms of the licence agreement.


☒ **I accept the terms of the licence agreement.***

Create Licence

Copyright © 1998-2020 Matrix Science Ltd. All Rights Reserved.

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Name = Human_EST Family = C:/inetpub/mascot/sequence/Human_EST_*_fasta
 Filename = Pathname =
 Status = Not in use
 State Time = Thu Jan 1 09:00:00 # searches = 0
 Mem mapped = NO Request to mem map = YES Request unmap = NO Mem locked = NO
 Number of threads = -1 Current = NO Type = Nucleic acid

Name = NIST_BSA_IonTrap Family = C:/inetpub/mascot/sequence/NIST_BSA_IonTrap/current/NIST_BSA_IonTrap_*.msp
 Filename = Pathname =
 Status = Not in use
 State Time = Thu Jan 1 09:00:00 # searches = 0
 Mem mapped = NO Request to mem map = YES Request unmap = NO Mem locked = NO
 Number of threads = -1 Current = NO Type = Nucleic acid

Mascot Server Licence Information

[Register a new product key](#) [View database status](#) [Reload this page](#)

Please include all the contents of this page when requested to provide this information to technical support.

Mascot Server version: 3.0.0
 Licence path: C:/inetpub/mascot/config/licdb

Licence(s) found:

Product Key	Start	End	Status
Active:			
Feature: Mascot Server - Core Functionality (v3.0)			OK
Feature: Mascot Server - CPU units (2)			
Feature: Mascot Server - Cores per CPU units (6)			
Company: MSKK			
User: Reimi SUZUKI			
Distributor: Matrix Science KK			
Inactive:			
Feature: Mascot Server - Core Functionality (v2.8)			OK
Feature: Mascot Server - CPU units (2)			
Feature: Mascot Server - Cores per CPU units (6)			
Company: MSKK			
User: Reimi SUZUKI			
Distributor: Matrix Science KK			
Feature: Mascot Server - Core Functionality (v2.7)			OK
Feature: Mascot Server - CPU units (2)			

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The screenshot shows the Mascot Server 3.0 web interface. The top navigation bar includes links: Home, Access Mascot Server, Help, Support, Training, Monthly newsletter, and About. The main content area displays the Mascot Server 3.0 logo and a description of the search engine. On the right side, there is a 'Quick Links' section with the following links: Log in, Search Log, Database Status, Configuration Editor, and Log out. The 'Configuration Editor' link is circled in red. Below this, there is a 'Database Security' section with links: Current session, Information, Change password, and Edit profile. A yellow box with the text 'Home の Configuration Editor -> Database Manager' is overlaid on the page. A large orange arrow points from the 'Configuration Editor' link in the Quick Links section to the 'Database Manager' link in the 'Mascot Configuration' table below.

Home の Configuration Editor -> Database Manager

Mascot Configuration

Amino Acids	Amino Acid Data
Modifications	Modification definitions
Symbols	Symbols used in chemical formulae
Linkers	Linker definitions
Enzymes	Enzyme definitions
Instruments	Fragmentation Rules
Quantitation	Quantitation Methods
Crosslinking	Crosslinking Methods
Configuration Options	Global Options in mascot.dat
Database Manager	Sequence databases, Parse Rules and automated downloads

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
Proxy settings

You may need to define proxy settings to enable Database Manager to access the Internet, depending on your local network and Internet connection.

Proxy settings define

- HTTP proxy str
- HTTPS proxy s
- FTP proxy str

HTTP test **succeeds**; can access http://http-connectivity-test.matrixscience.com/db_manager.txt from the Mascot server machine.

HTTPS test **succeeds**; can access https://https-connectivity-test.matrixscience.com/db_manager.txt from the Mascot server machine.

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Proxy settings

mascot.dat settings

The recommended proxy type is Auto, or failing that, Environment, Registry, Specify, or WPAD. Auto means Database Manager will try to find proxy settings from the system registry.

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https://www.matrixscience.com/help/seq_db_setup.html

https://www.matrixscience.com/help/seq_db_setup_library.html

cRAP	Enable	NIST_Drosophila_IonTrap	Enable
Environmental_EST	Enable	NIST_E.coli_IonTrap	Enable
Fungi_EST	Enable	NIST_HSA_IonTrap	Enable
Human_EST	Enable	NIST_Human_HCD	Enable
Invertebrates_EST	Enable	NIST_Human_HCD_2_good	Enable
Mammals_EST	Enable	NIST_Human_HCD_3_semitryp	Enable
Mus_EST	Enable	NIST_Human_HCD_iTRAQ_1	Enable
NCBIprot	Enable	NIST_Human_HCD_iTRAQ_2	Enable
neXtProt	Enable	NIST_Human_HCD_iTRAQ_Phospho	Enable
Plants_EST	Enable	NIST_Human_IonTrap	Enable
Prokaryotes_EST	Enable	NIST_Mouse_HCD	Enable
Rodents_EST	Enable	NIST_Mouse_HCD_iTRAQ	Enable
SwissProt_AC	Enable	NIST_Mouse_HCD_iTRAQ_Phospho	Enable
SwissProt_ID	Enable	NIST_Mouse_IonTrap	Enable
Trembl_AC	Enable	NIST_Rat_IonTrap	Enable
Trembl_ID	Enable	NIST_Rat_QToF	Enable
UniRef100	Enable	NIST_S.cerevesiae_IonTrap	Enable
UP186698_X_laevis	Enable	NIST_S.cerevesiae_QToF	Enable
UP1940_C_elegans	Enable	NIST_S.pombe_IonTrap	Enable
UP2195_D_discoideum	Enable	NIST_UPS1_IonTrap	Enable
UP2311_S_cerevisiae	Enable	PRIDE_Arabidopsis	Enable
UP2485_S_pombe	Enable	PRIDE_B.burgdorferi	Enable
UP2494_R_norvegicus	Enable	PRIDE_B.cereus	Enable
UP437_D_rerio	Enable	PRIDE_B.subtilis	Enable
UP5226_T_rubripes	Enable	PRIDE_C.elegans	Enable
UP5640_H_sapiens	Enable	PRIDE_Contaminants	Enable
UP589_M_musculus	Enable	PRIDE_Drosophila	Enable
UP59680_O_sativa	Enable	PRIDE_E.coli	Enable
UP625_E_coli_K12	Enable	PRIDE_G.sulfurreducens	Enable
UP6548_A_thaliana	Enable	PRIDE_Human	Enable
UP6906_C_reinhardtii	Enable	PRIDE_Mouse	Enable
UP7305_Z_mays	Enable	PRIDE_Rat	Enable
UP803_D_melanogaster	Enable		
UP808_M_pneumoniae	Enable		
UP9136_B_taurus	Enable		

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The screenshots show the following steps:

- Step 1:** A list of predefined databases with 'Enable' buttons. The 'NIST_Human_IonTrap' button is highlighted with a red box.
- Step 2:** The 'Enable predefined definition' screen. The 'Next' button is highlighted with a red box.
- Step 3:** The 'Enable' screen for the 'NIST_Human_IonTrap' database. The 'Create' button is highlighted with a red box.
- Step 4:** The 'Summary' screen showing the configuration details for the 'NIST_Human_IonTrap' database. The 'Create' button is highlighted with a red box.

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<対象生物種リスト:DB名XXXXに該当> Environmental , Fungi, Invertebrates, Mammals, Mus, Plants, Prokaryotes, Rodents, Vertebrates

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The screenshot shows the Mascot Database Manager interface. On the left, the 'Running tasks (1)' link is highlighted. The main panel displays 'Database Manager tasks' with a table of tasks. A yellow box highlights the 'Task log' section, showing a successful update for 'PRIDE_Contaminants'.

Task	Status	Message
Download files for SwissProt	Running	[3.5%] Downloading /human_hcd_select
Download files for PRIDE_Contaminants	Completed	(success) 'PRIDE_Contaminants' successfully updated.

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Name	= PRIDE Contaminants	Family	= C:/inetpub/mascot/sequence/PRIDE Cont
Filename	= PRIDE Contaminants_20160906.msg	Pathname	= C:/inetpub/mascot/sequence
Status	= not in use	Statistics	
State Time = Fri Oct 4 11:49:45 # searches = 0			
Mem mapped = NO Request to mem map = NO Request unmap = YES Mem locked = NO			
Number of threads = 1 Current = NO Type = Spectral library			
Name	= PRIDE Contaminants	Family	= C:/inetpub/mascot/sequence/PRIDE Cont
Filename	= PRIDE Contaminants_20160908.msg	Pathname	= C:/inetpub/mascot/sequence
Status	= In use	Statistics	Compression warnings Recompress file
State Time = Fri Oct 4 11:49:45 # searches = 0			
Mem mapped = NO Request to mem map = NO Request unmap = NO Mem locked = NO			
Number of threads = 1 Current = YES Type = Spectral library			
Name	= SwissProt	Family	= C:/inetpub/mascot/sequence/SwissProt/current/SwissProt *.fasta
Filename	= SwissProt_2024_04.fasta	Pathname	= C:/inetpub/mascot/sequence/SwissProt/current/SwissProt_2024_04.fasta
Status	= In use	Statistics	Unidentified taxonomy Recompress file
State Time = Fri Oct 4 11:39:45 # searches = 0			
Mem mapped = YES Request to mem map = YES Request unmap = NO Mem locked = NO			
Number of threads = -1 Current = YES Type = Amino acid			
Name	= SwissProt	Family	= C:/inetpub/mascot/sequence/SwissProt/current/SwissProt *.fasta
Filename	= SwissProt_2024_05.fasta	Pathname	= C:/inetpub/mascot/sequence/SwissProt/current/SwissProt_2024_05.fasta
Status	= Creating compressed files 1% complete		
State Time = Fri Oct 4 12:29:19 # searches = 0			
Mem mapped = NO Request to mem map = YES Request unmap = NO Mem locked = NO			
Number of threads = -1 Current = NO Type = Amino acid			

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